



**ALMOSTAKBAL**  
For Dental Supplies

المستقبل لمستلزمات طب الأسنان



**JOIN  
CHAMP**

# ZC-A200 NEW



**ALMOSTAKBAL**  
For Dental Supplies  
المستقبل لمستلزمات طب الأسنان



# JOIN CHAMP ZC-A300



**ALMOSTAKBAL**  
For Dental Supplies  
المستقبل لمستلزمات طب الأسنان



JOIN  
CHAMP

# ZC-S400 NEW



ALMOSTAKBAL  
For Dental Supplies  
المستقبل لمستلزمات طب الأسنان



# JOIN CHAMP ZC-A500



**ALMOSTAKBAL**  
For Dental Supplies  
المستقبل لمستلزمات طب الأسنان



# JOIN CHAMP ZC-S500



**ALMOSTAKBAL**  
For Dental Supplies  
المستقبل لمستلزمات طب الأسنان





**Configuration**

PU cushion(0.9mm)

Double arm rest

LED reflective lamp

6 in 1 foot pedal

Assistant control system USA electric motor

90 ° Rotation cabinet box, convenient four-hand operation

German water & air tubes

Cup Holders and Trays



Yadens

# YD-A5(New)



ALMOSTAKBAL  
For Dental Supplies

المستقبل لمستلزمات طب الأسنان



- Configuration
- 6 beads led lamp
- Double arm rest
- 24V film view
- Instrument silicon pad
- PU cushion(0.9mm)
- 90° Rotation cabinet box,convenient four-hand operation
- USA electric motor
- German water& air tubes
- Assistant control system
- 3 in 1 tissue box
- 6 in 1 foot pedal



## YD-D1



### Configuration:

- Streamlined chair design
- PU Cushion 1.2mm
- Double arm rest
- Stainless foot pedal tubes
- USA electric motor
- German water&air tube
- Instrument silicon pad and assistant silicon pad
- 6 in 1 foot pedal
- 024 doctor stool
- Assistant control system with silicon pad.
- 6 beads direct Led lamp
- Stainless steel foot pedal outer tube
- Import handpiece inner tubes and stainless connection
- Copper water supply pipe





1. 6 in 1 foot pedal
2. 041 luxury doctor stool
3. Assistant control system with silicon pad.
4. Casting Aluminum luxury led lamp
5. 2.2N Luxury lamp arm
6. Stainless steel foot pedal outer tube
7. Double water bottle
8. Linkage chair frame
9. Detachable casting aluminum backrest
10. 4.3" LCD display screen
11. Import handpiece inner tubes and stainless connection
12. Mirco-Fiber leather
13. Aluminum headrest
14. 9 memories programs
15. Silicon strong & weak suction
16. Package with accessories:  
3-in-1 Tissue Bucket





**Configuration:**

1. Casting aluminum construction for assistant, cabinet box, instrument base
2. PBT injection molding filter
3. 90° Rotation cabinet box
4. pillow for head and waist
5. Double water bottle
6. 4.3" main control LCD screen with 9 memories programs
7. Stainless steel handpiece connector
8. Micro-fiber leather
9. G8 led lamp
10. 072 Dentist stool and 026 Nurse Stool
11. GL-G8 8 beads luxury led lamp
12. Double Armrest
13. 6 in 1 foot pedal
14. Square lamp arm
15. Import Noiseless motor
16. Import handpiece inner tubes
17. Stainless foot pedal tubes





### Configuration:

1. Casting aluminum construction for assistant, cabinet box, instrument base
2. PBT injection molding filter
3. 90° Rotation cabinet box
4. Pillow for head and waist
5. Double water bottle
6. 7" main control LCD screen and 4.3" assistant control LCD screen.
7. 9 memories programs
8. Stainless steel handpiece connector
8. Real leather cushion
9. G8 LED lamp
10. 072 Dentist stool and 026 Nurse Stool
11. GL-G8 8 beads luxury LED lamp
12. Double Armrest
13. 6 in 1 foot pedal
14. Square lamp arm
15. Import Noiseless motor
16. Import handpiece inner tubes
17. Stainless foot pedal tubes

**AutoClave**

**SKY**



**Class B  
18 Liter**



**Class B  
23 Liter**



**Class N  
18 Liter  
Digital**



**Class N  
18 Liter  
5 Programs**



**AutoClave**

**Accessories**



**Sealing Machine  
Metal**



**Sealing Machine  
Fiber**



**Water Distiller**

## Technical Datasheet Class B Autoclave



### Commodity Parameter:

Voltage: AC 220V/110v ±10%

Frequency: 50/60Hz

Wire plug: GB 3 cores/EN 3 cores/ American

Max Hp: 2000VA

Tank Capacity: 18/ 23L

Class No: European Class B standard

Sterilizing temperature: 134°C/121°C, 136-138 Optional

Inside Dimension: Φ249×450mm

Water tank : SS304 wide open on top easy for cleaning

Display board: large LED Screen

## Technical Datasheet Class B Autoclave

### Technical Parameters:

Sterilization level : European Class B standard

Sterilizing temperature: standard 134°C/121°C.( 136-138 °C adjustable))

Special sterilization: Killing active HV,HBV,BSE and spore,

Dry procedure: Strong vacuum drying ,the residual humidity<0.2%

Display: Widescreen LED showing time, temp and pressure and real time operation procedure.

Test: Bowie & Dick Steam penetration test, vacuum test , Helix test

Safety: EU standard Safety valve,

Lock Security: Double door locking security, convenient and users friendly..

Sterilizing record: Optional USB and printer.

Water supply system: 3L distiller water tank,3L waste water tank.

Instrument container : 5 levels with 3 mesh trays.

Both Clean water shortage alarming and waste water overloaded warning system.



### Product Characters:

1. Based on the new class B of international prEN13060, have 3 pulsating vacuum, Vacuum up to -0.95Bar (-0.09Mpa). Suitable to sterilizing wrap, unwrap, solid, porous, hollow device.
2. Use imported & advanced 16-bit microprocessor, easy to operate. Most suitable for the Dept of Stomatology , Ophthalmology, Surgery and Lab.
3. With B& D test and vacuum test procedures. To test the penetrability.
4. Chamber punched by Stainless Steel (#304, thickness: 2.5mm), together with The jet type of steam generator, self-cleaning ensure the efficient sterilizing.
5. With alarm system for waste water tank, avoid the waste water to enter into the sterilizing circulation, make the sterilizing completely.
6. Precise LCD digital display and advanced self-test system ensure you get the working data freely.
7. Double lock door system for safety during operating.
8. With built-in printer .USB optional. The disinfection process data and date automatically print.

# Portable Dental X-Ray

RV-210RG



Adjustable KV & Current

Ultra-low Radiation Leakage

**SKY** Ray

## Features



- Single-person operation:**  
Designed for easy operation by a single user.
- Portable design:**  
Compact and lightweight for convenient use in various settings.
- Adjustable imaging angles:**  
Provides flexibility to adjust imaging angles for optimal results.
- Efficient image capture:**  
Enables quick and efficient capture of dental X-ray images.
- User-friendly interface:**  
Intuitive controls make it easy to operate.
- Versatile applications:**  
Suitable for a wide range of dental imaging needs.
- Safety features:**  
Incorporates safety measures to minimize radiation exposure.
- Reliable performance:**  
Delivers consistent and reliable imaging results.

## Specifications

Item	Specification
Tube Voltage	60kV to 70kV adjustable
Tube Current	2.5 mA, 3.5 mA adjustable
KV Ripple Frequency	< 1%
Inverting Frequency	> 100kHz
Screen	Capacitive touch screen
Step of Exposure Time	R20 Coefficient
Exposure Time	0.01 – 1 Second

# HIGH-FREQUENCY PORTABLE DENTAL X-RAY MACHINE R-RAY I

**Accurate imaging:** X-ray machines can produce high-quality images that allow medical professionals to diagnose and treat a wide range of conditions.  
**Non-invasive:** Unlike other medical procedures, X-rays are non-invasive and do not require any incisions or injections.  
**Quick results:** X-rays can provide results quickly, allowing medical professionals to make timely decisions about treatment.



Pure imported Canon high-frequency ray tube



Stable high-frequency inverter output technology



Ultra high capacity customized lithium batteries



RVG digital mode one click conversion



Automatic error reporting for routine POST faults



Humanized color LCD touch setting interface



Built in temperature control, overvoltage, undervoltage protection

## Dental Imaging System



## Dental Digital Sensor



Model  
Sensor dimension  
Imaging Area  
Pixel Matrix  
Sensor technology  
Scintillator  
Length of cable  
Pixel pitch  
Theoretical resolution  
Grayscale  
X-ray trigger mode  
Ingress protection of sensor  
Additional function  
OS

**Sensor 1.5**  
29.2x42.5x4.5  
24 x 33 mm  
(1.98M )1200 x 1650  
APS CMOS  
Csl (TI) without FOP  
2.7m  
20µm  
25 lp/mm  
65,535 gray level(16 bit)  
AED  
IP68 for Soak disinfection  
support TWAIN protocol  
Windows 7/10/11, 32/64 bits

## Dental Compressor

## Dental Noiseless Oilless Air Compressor



- Oil Free
- 40 Litre



- Oil Free
- 50 Litre



- Oil Free
- 40 Litre

## Bleaching Machine



**ATLAS**  
14 Lamps

Details:

input power: 100-240v, 50/60Hz

power: 4-46w

wave length: 395-625mm

Density: 200-300mw/cm<sup>2</sup>

packing weight: 2.5 Kg

adjustable timer: 130 min

light resource: 14 pcs high power LED



## Multi Media



**intra Oral Camera  
With Monitor Arm**



**HD- Intraoral Camera**



**Micromirror  
Camera+Display**



## Light Cure

**SIVEN**  
FOR FUTURE DENTISTRY  
**LED Light Cure  
One Second**



**iMax  
One Second**



**LED B  
Light Cure**





**VG-MP** 迷你按压手机  
Mini head push button handpiece

**VG-MPL** 迷你风光按压手机  
Mini head LED push button handpiece



**VG-SP2** MAX-2标准按压手机  
Standard push button handpiece



**VG-TP3-3** MAX-3大头按压手机  
MAX-3 Super torque push button handpiece



**VG-TP** MAX-1大头按压手机  
MAX-1 Super torque handpiece



**VG-TP-2** 大头按压手机  
Super torque handpiece (type 2)



**VG-TP3-1** 大头按压手机(三点喷水)  
Super torque handpiece with triple water spray (type 1)



**VG-TP3-2** 大头按压手机(三点喷水)  
Super torque handpiece with triple water spray (type 2)

### 技术参数 / Technical

取针方式/ Chuck type	Push button chuck
转速/Rotation speed	≥300,000 rpm
工作气压/Air pressure	0.2-0.25Mpa(2-2.5kgf/cm <sup>2</sup> )
适用车针/Bur applicable	1.59mm-1.6mm
噪音/Noise	≤68dB



VG-5TPL

五孔大头风光手机  
5 LED super torque LED handpiece



VG-TPL

大头风光手机(款式1)  
Super torque LED handpiece (type 1)



VG-TPL2

大头风光手机(款式2)  
Super torque LED handpiece (type 2)

### 技术参数 / Technical

取针方式/ Chuck type	Push button chuck
转速/Rotation speed	300,000-350,000 rpm
工作气压/Air pressure	0.2-0.25Mpa(2-2.5kgf/cm <sup>2</sup> )
适用车针/Bur applicable	0.1-0.25Mpa
噪音/Noise	≤68dB



VG-45SP

45度角高速手机  
45° High speed handpiece



VG-45SPL

45度角高速手机  
45° LED handpiece

### 技术参数 / Technical

气压/Air pressure	0.2-0.25Mpa(2-2.5kgf/cm <sup>2</sup> )
水压/Water pressure	0.1Mpa-0.25Mpa
转速/Rotation speed	≥300,000rpm
适用车针/Bur applicable	1.59mm-1.6mm
噪音/Noise	≤68dB

VIGOR | 01



VG-SH01



VG-SH02



VG-CA01



VG-CA02



VG-UCA02



VG-AM01



VG-AM02

### E型低速手机套装(外置水道) / E type low speed handpiece (External cooling system)

气压/Air pressure	≤245-395Kpa
转速/Rotation speed	18,000-22,000 rpm
适用车针/Bur applicable	Ø2.35mm
噪音/Noise	≤70dB

**TorNa**<sup>TM</sup>  
Endo



**ALMOSTAKBAL**  
For Dental Supplies

المستقبل لمستلزمات طب الأسنان

# DENTAL DAM

- ✓ Better Colour Contrast To Enhance Tooth Visibility
- ✓ Powder Free
- ✓ High Stretchability
- ✓ High Tear Resistance
- ✓ Endodontic treatment.
- ✓ Protect the treated area from being exposed to oral bacteria and salivary contamination.
- ✓ Prevent the patient's tongue and gums from irritation during treatment.
- ✓ Prevent any debris and or dental instruments from falling to the patient's throat.

# T-ONE

- Small and cute in design and portable
- High accuracy
- Automatic calibration
- High resolution and big screen
- Adjustable viewing angle

Portable, High accuracy, Automatic calibration, Large battery capacity



## ENDO Motor Dolphin



Product is based on ergonomic design with a lightweight endomotor  
 There are 6 kinds of procedures (Memory Set) for selection. It can input the optimum data in advance, and also can modify it accordingly to the working conditions in use. Dental Contra-angle handpiece and lamp holder can withstand 135C high temperature and high pressure sterilisation  
 In-built high capacity lithium battery, so it can reduce battery change times and reducing costs of the machine.





**ALMOSTAKBAL**  
For Dental Supplies  
المستقبل للمستلزمات طب الأسنان



**Premium LIGHT BODY**

- Secondary precision impression material for abutment tooth
- Inlay, Onlay, Crown and Bridge



**HanFil Flow**

- Excellent flow
- Low shrinkage
- High strength and polishability
- Natural esthetics shade



**Han Luting Cement**

- Dual cured system
- Fast and safe applications
- Self Etching and Self adhesive
- Easy clean up of excess



**Premium PUTTY**

- Primary impression for tray bases
- Impression for two-step techniques



**HanBase Liner**

- Primary impression for tray bases
- Impression for two-step techniques



**Han Temp Cement NE**

- Suitable bond strength
- Easy removal of the cement
- Formulated free of Eugenol
- Prevent sensitivity to cold and heat



**Han Temp Cement EU**

- Effect for dental pulp with Eugenol
- Suitable bond strength
- Easy removal of the cement
- Prevent leakage of oral fluids



**HanTemp**

- High elasticity and strength
- Easy to remove
- Tight sealing effect
- Easy application



**HanFil Composite**

- High strength and polishability
- Excellent handling
- Low shrinkage
- Excellent matching with natural tooth



**HanFil Composite Kit**

- HanFil 4g x 5 Syringes
- HanEtch 3ml x 2 Syringes
- HanBond 5ml x 1 Bottle
- Accessory ; Micro brush



**Core build-up**

- Easy application
- Dual cured system
- High compressive strength
- High flexural strength



**HanTemp Crown**

- Easy handling and nice injection
- Low shrinkage
- High strength and abrasion resistance
- Protect a dental pulp



**HanEtch-37**

- HanFil 4g x 5 Syringes
- HanEtch 3ml x 2 Syringes
- HanBond 5ml x 1 Bottle
- Accessory : Micro brush



**Han Bond 7**

- No need etchant
- Light curing dentin adhesive
- Easy and fast application
- No stickiness and good flow
- 7th generation



FOR BETTER  
ENDO

**NIC<sup>®</sup>**



# ENDODONTIC INSTRUMENTS



**ALMOSTAKBAL**  
For Dental Supplies

المستقبل لمستلزمات طب الأسنان

# EZ-Pass


Launch of NIC EZ-PASS®, the established crown-down NiTi system. EZ-PASS® is applicable to 80% regular cases.



11mm 21/25/31mm


## CHARACTERISTIC

- MaxTech 6.0 NiTi alloy
- Fixed tapered for easy obturation
- Convex triangle cross-section exhibits excellent centricity
- Progressive screw pitches increase debris removal and less screw-in effect



## OPERATING STEPS

EZ-PASS						
Model	19mm	21mm	25mm	25mm	31mm	31mm
● 20/.07	EZP 031 2007	-	-	-	-	-
○ 15/.04	-	EZP 021 1504	EZP 025 1504	EZP 031 1504	-	-
● 20/.04	-	EZP 021 2004	EZP 025 2004	EZP 031 2004	-	-
● 20/.06	-	EZP 021 2006	EZP 025 2006	EZP 031 2006	-	-
● 25/.04	-	EZP 021 2504	EZP 025 2504	EZP 031 2504	-	-
● 25/.06	-	EZP 021 2506	EZP 025 2506	EZP 031 2506	-	-
Sorted	-	EZP 021 A001	EZP 025 A001	EZP 031 A001	-	-



## USAGE STEPS

- Prepare straightline access to canal orifice
- Explore the canal with small size hand files(e.g. SS K-file #10)
- Open 2/3 upper coronal part with opener P1 20/.07
- Shape and finish the canal with P2 15/.04, P3 20/.04, P4 20/.06, P5 25/06 respectively to reach W/L passively by brushing
- EDTA, Irrigation is recommended to be used for better result

Model	19mm	21mm	25mm	31mm
● 20/.07	EZP 031 2007	-	-	-
○ 15/.04	-	EZP 021 1504	EZP 025 1504	EZP 031 1504
● 20/.04	-	EZP 021 2004	EZP 025 2004	EZP 031 2004
● 20/.06	-	EZP 021 2006	EZP 025 2006	EZP 031 2006
● 25/.04	-	EZP 021 2504	EZP 025 2504	EZP 031 2504
● 25/.06	-	EZP 021 2506	EZP 025 2506	EZP 031 2506
Sorted	-	EZP 021 A001	EZP 025 A001	EZP 031 A001



STERILE R

# EZ-Large

EZ-LARGE is designed for large root canal preparation. Doctors can choose the appropriate model based on the clinical root canal cases.



11mm 21/25/31mm

EZ-LARGE						
Model	21mm	25mm	31mm			
L1	EZL 021 3004	EZL 025 3004	EZL 031 3004			
L2	EZL 021 3006	EZL 025 3006	EZL 031 3006			
L3	EZL 021 3504	EZL 025 3504	EZL 031 3504			
L4	EZL 021 3506	EZL 025 3506	EZL 031 3506			
L5	EZL 021 4004	EZL 025 4004	EZL 031 4004			
L6	EZL 021 5004	EZL 025 5004	EZL 031 5004			



Model	21mm	25mm	31mm
L1	EZL 021 3004	EZL 025 3004	EZL 031 3004
L2	EZL 021 3006	EZL 025 3006	EZL 031 3006
L3	EZL 021 3504	EZL 025 3504	EZL 031 3504
L4	EZL 021 3506	EZL 025 3506	EZL 031 3506
L5	EZL 021 4004	EZL 025 4004	EZL 031 4004
L6	EZL 021 5004	EZL 025 5004	EZL 031 5004



nicmedical.com

# GF File®

Patented Technology  
Gradient flexibility  
Archimedes' spiral movement



# GF File



11mm 21/25/31mm

## CHARACTERISTIC

- Gradient flexibility is produced with patented technology, presenting different flexibility.
- Progressive flexibility transition from the upper to the apical.
- Depth marking area helps the file to be well centered in the root canal and effectively conduct axial force source.
- Respect natural root canal anatomy, get accurate feedback from file tip with true expression of root canal morphology.
- Sequential shaping reduce the difficulty of preparing small and curved root canals.
- One last standardized finishing file Z3 25/.04 for intend to taper confirmation and precise match to GP.

STERILE R



## PARAMETERS

GF FILE					
Model	19mm	21mm	25mm	25mm	31mm
● SK	GFF 019 1504	-	-	-	-
○ T	-	GFF 021 1702	GFF 025 1702	GFF 031 1702	-
● Z1	-	GFF 021 2004	GFF 025 2004	GFF 031 2004	-
● Z2	-	GFF 021 2006	GFF 025 2006	GFF 031 2006	-
● Z3	-	GFF 021 2506	GFF 025 2506	GFF 031 2506	-
SK-Z3	-	GFF 021 A001	GFF 025 A001	GFF 031 A001	-

Model	19mm	21mm	25mm	31mm
● SK	GFF 019 1504	-	-	-
○ T	-	GFF 021 1702	GFF 025 1702	GFF 031 1702
● Z1	-	GFF 021 2004	GFF 025 2004	GFF 031 2004
● Z2	-	GFF 021 2006	GFF 025 2006	GFF 031 2006
● Z3	-	GFF 021 2506	GFF 025 2506	GFF 031 2506
SK-Z3	-	GFF 021 A001	GFF 025 A001	GFF 031 A001



nicmedical.com

# K-Flex



Model	21mm	25mm	28mm	31mm
○ 15	TFK 021 015	TFK 025 015	TFK 028 015	TFK 031 015
● 20	TFK 021 020	TFK 025 020	TFK 028 020	TFK 031 020
● 25	TFK 021 025	TFK 025 025	TFK 028 025	TFK 031 025
● 30	TFK 021 030	TFK 025 030	TFK 028 030	TFK 031 030
● 35	TFK 021 035	TFK 025 035	TFK 028 035	TFK 031 035
● 40	TFK 021 040	TFK 025 040	TFK 028 040	TFK 031 040
15-40	TFK 021 A01	TFK 025 A01	TFK 028 A01	TFK 031 0A01

# H-File



NIC H-Files are used for probing, permeabilization or extraction of debris. Their use is essentially for widening after passage of the K-file of the same number.



## CHARACTERISTIC

- Available in Niti Alloy and Stainless Steel material
- Round cross section
- Autoclave

STERILE R




## S.S. Hand Use H-File

Size	21mm	25mm	28mm	31mm
6	HSH 021 006	HSH 025 006	HSH 028 006	HSH 031 006
8	HSH 021 008	HSH 025 008	HSH 028 008	HSH 031 008
10	HSH 021 010	HSH 025 010	HSH 028 010	HSH 031 010
15	HSH 021 015	HSH 025 015	HSH 028 015	HSH 031 015
20	HSH 021 020	HSH 025 020	HSH 028 020	HSH 031 020
25	HSH 021 025	HSH 025 025	HSH 028 025	HSH 031 025
30	HSH 021 030	HSH 025 030	HSH 028 030	HSH 031 030
35	HSH 021 035	HSH 025 035	HSH 028 035	HSH 031 035
40	HSH 021 040	HSH 025 040	HSH 028 040	HSH 031 040
45	HSH 021 045	HSH 025 045	HSH 028 045	HSH 031 045
50	HSH 021 050	HSH 025 050	HSH 028 050	HSH 031 050
55	HSH 021 055	HSH 025 055	HSH 028 055	HSH 031 055
60	HSH 021 060	HSH 025 060	HSH 028 060	HSH 031 060
70	HSH 021 070	HSH 025 070	HSH 028 070	HSH 031 070
80	HSH 021 080	HSH 025 080	HSH 028 080	HSH 031 080
15-40	HSH 021 A01	HSH 025 A01	HSH 028 A01	HSH 031 A01
45-80	HSH 021 A02	HSH 025 A02	HSH 028 A02	HSH 031 A02

## NiTi Hand Use H-File

Size	21mm	25mm	28mm	31mm
6	HNH 021 006	HNH 025 006	HNH 028 006	HNH 031 006
8	HNH 021 008	HNH 025 008	HNH 028 008	HNH 031 008
10	HNH 021 010	HNH 025 010	HNH 028 010	HNH 031 010
15	HNH 021 015	HNH 025 015	HNH 028 015	HNH 031 015
20	HNH 021 020	HNH 025 020	HNH 028 020	HNH 031 020
25	HNH 021 025	HNH 025 025	HNH 028 025	HNH 031 025
30	HNH 021 030	HNH 025 030	HNH 028 030	HNH 031 030
35	HNH 021 035	HNH 025 035	HNH 028 035	HNH 031 035
40	HNH 021 040	HNH 025 040	HNH 028 040	HNH 031 040
45	HNH 021 045	HNH 025 045	HNH 028 045	HNH 031 045
50	HNH 021 050	HNH 025 050	HNH 028 050	HNH 031 050
55	HNH 021 055	HNH 025 055	HNH 028 055	HNH 031 055
60	HNH 021 060	HNH 025 060	HNH 028 060	HNH 031 060
70	HNH 021 070	HNH 025 070	HNH 028 070	HNH 031 070
80	HNH 021 080	HNH 025 080	HNH 028 080	HNH 031 080
15-40	HNH 021 A01	HNH 025 A01	HNH 028 A01	HNH 031 A01
45-80	HNH 021 A02	HNH 025 A02	HNH 028 A02	HNH 031 A02



# Gates Drill



Gates drill used for expanding the canal entrance or preparation for middle-upper part of the canals



## CHARACTERISTIC

- Medical Stainless Steel
- Sizes: #1, #2, #3, #4, #5, #6
- Length: 28mm - 32mm
- Recommended speed: 800-1200rpm
- Corresponding sizing rings



Size	Rings	28mm	32mm
1 / 050	1	ESG 028 001	ESG 032 001
2 / 070	2	ESG 028 002	ESG 032 002
3 / 090	3	ESG 028 003	ESG 032 003
4 / 110	4	ESG 028 004	ESG 032 004
5 / 130	5	ESG 028 005	ESG 032 005
6 / 150	6	ESG 028 006	ESG 032 006
1-6	-	ESG 028 A01	ESG 032 A01



# SIVEN®

FOR FUTURE DENTISTRY

## MORE THAN ENDO

QUALITY SURE, CUSTOMIZED SERVICE, ONE-STOP SOLUTION

Dental Root Canal Instruments



SV-Kids Files



SV-Flexi



SV-ONE



SV-TAPER




SV-SHARPER



D Files



C Files



Path Files



## Instruments



**Orthodontic Pliers**



**Extraction Forceps**



**Rubber Dam Instruments**

## ROOT ELEVATOR



## Orthodontic



Power Chain



Titanium Micro Implant

## Elastics



## Brackets & Wires



Wires



Face Mask





## Disposable



**Sectional Matrices**



**Brushes**



**Retraction Cord**



**Irrigation Needle**



**Lubricant**



**Sutures**



**Composite Polishing Kit**

# Disposable Burs & Stones



**Matrice Bands**



**Articulating Paper**



**X-Ray Film Holders**

**Irrigation Tips**

**Brushes**

**Applicators**



## مركز صيانة المستقبل

- نتشرف أن نعلن لعملائنا الكرام عن قسم الصيانة الخاص بمؤسسة المستقبل للإستيراد و التصدير .

- لدينا فريق صيانة علي أعلى مستوى يقدم للعملاء خدمة ما بعد البيع بأسلوب حديث حيث أن فريق العمل مدرب و لديه شهادات معتمدة .

- لدينا قطع غيار وحدات الأسنان و جميع الأجهزة التي تستوردها الشركة .

- لدينا إمكانية تجديد عيادات الأسنان و تبديل الوحدات القديمة بوحدات حديثة .

- تضمن الشركة لعملائها الكرام السرعة و الجدية في صيانة أجهزتها في وقت سريع جدا .

الصيانة : 01050698355 01050698311



We are building the future

SCAN ME



العنوان : مصر - الدقهلية - المنصورة  
شارع الأشقر أمام مستشفى الكلى

الإدارة : 01094961076 01009041812  
01000322069

المبيعات : 01050119944 01050119933  
01050698344

E-MAIL : waleedsaad9300@gmail.com