

# MEDiDENT

Innovative Dental Solutions


*Unseen Power Beneath The Surface*

Tapered • Trusted • Titanium

# Surface Treatment Ultra SLA (Active)

**Highness**  
**HSN**  
Hydrophilic Nano Texture

- Surface treatment of the noblesse line is performed on 2 steps ,first SLA , then Alkaline Cleansing ,providing a hydrophilic nano and microscale texture improving osseointegration.
- 11° Morse taper connection provides a stable and predictable abutment implant connection.
- Highness HS1 Design provides a progressive thread design for optimum stress distribution.
- Highness HSVII design provides aggressive thread design for cases with soft bone for optimum primary stability.



## Implant Platform

Platform switching between implant and abutment improve soft tissue healing and preserves crestal bone loss



## Ultra SLA Surface

A combination of SLA followed by Alkaline cleansing allows for a hydrophilic micro and Nano surface texture improving osseointegration and BIC

## Wide Cutting Edge

self-cutting edge for enhanced primary stability & improved anti rotation



## 2.1Hex, 2.5Hex

Excellent compatibility with domestic and foreign products

## Square Thread

Wide square threads, Allows for compressive stress distribution and preserves crestal bone

## Active Threads

Deep , sharp threads to maximize initial stability even for immediate loading opportunities and long-term success contact

## Implant Body

Tapered implant body

## Apex Thread

Blunt Apex for safe implant insertion



# Runyes Family

## Runyes 3DS V5



### Expansive Vision & Depth

16×14mm field of view, 25mm depth for efficient, precise scanning.

### High Accuracy

10µm single tooth, 20µm full arch.

### Faster Scanning

Full scan in 60s with AI algorithms.

### Lightweight Design

160g ergonomic handpiece for effortless use.

### Gesture Control

View models with simple hand movements.

## Runyes 3DS V6

### Wireless Dual-Band WiFi

5+2.4 GHz for fast, stable data transmission.

### Long Battery Life

Dual batteries, up to 2 hours continuous scanning.

### Motion Gesture Sensor

View models with simple hand movements.

### Emergency Mode

Switch to wired connection anytime.

### Faster Scanning

AI algorithms, full scan in as fast as 60 seconds.

### High Accuracy

10 µm (single tooth), 20 µm (full arch).

### Lightweight

Only 230g including batteries, easy for long use.

### Wide & Deep View

16×14 mm FOV, 25 mm depth for efficient, precise scans.



## Photogrammetry Camera

### High Accuracy

20S Fast Half-Arch Impression

More Comfortable

High Stability

Simple Impression Workflow

Data Openness

No Need to Calibrate for Each Use Simple Operation

1.3kg Lightweight Design

Light Design for The Bracket

Making It Easy to Pick Up and Put Down

## Color Shade Detector

### Greater Accuracy

7 optical fiber heads for full tooth surface measurement, reducing errors.

### True Color Accuracy

Multiple algorithms, determines color in 1 second with one click.

### External Factor Resistant

LED white light eliminates effects of lighting and environment.

### Simple Interface

Clear, easy-to-read two-color display.

### Lightweight Design

216g handle for comfortable, long-term use.

### Ultra-long Battery

Smart standby, over 5000 measurements per charge.



Highness Implant Fixture Surface Feature

# Highness HSN

Hydrophilic Nano Fixture



## Contact us

- +2 010 120 014 99
- medidentegypt@gmail.com
- Head Office : Canan Capital Business centre,  
Office 734 , Katamya - Cairo



SCAN ME

Follow us Now



highness

Runyes  
FOR ALL DENTAL

Noblesse

aquapick

UBGEN

SANI