

# CASTELLINI

PASSION FOR DENTISTRY

Since 1935



PUMA ELI R




SKEMA 6



Made in Italy



**B CLASSIC**



**MILLSEAL+ MANUAL**



**TETHYS H10 PLUS**

# Corpus Vac®

Vacuum System and Equipment



DIRECT - SEMI WET






MAXI SUCTION



MIDI SUCTION

# clmp



# Kuraray



# Bossklein®

## Complete hygiene solutions



### DAILYasp

Aspirator Cleaner & Disinfectant

- Concentrated solution
- EN14476 virus tested
- 15 min contact time\*
- 20ml/1L dosage (2%)\*



### ULTRAbio

Ultrasonic Cleaner

- Concentrated solution
- Enzymatic
- 15 min contact time\*
- 30ml/1L dosage (3%)



### SDactiv ZERO

Surface Disinfectant

- Ready to use
- Alcohol free
- EN14476 tested



### V-WIPE ZERO

Surface Disinfectant Wipe

- Alcohol free
- EN14476 tested
- 60 second contact time
- Tubs and Refill Packs
- Three wipe types



### SDactiv

Surface Disinfectant

- Alcohol (63%)
- EN14476 tested
- 60 second contact time



